

# CRITIKON™ Blood Pressure Cuff Adapters and Connectors

Total solutions for non-invasive blood pressure



# GE Healthcare's connection system is simple and easy to implement.

GE Healthcare provides cost-effective standardization and inventory management solutions. By standardizing blood pressure (BP) connection systems, you create a universal connection to allow the same type of BP cuff connection to work on a wide variety of non-invasive blood pressure devices. By consolidating all of your BP cuff purchases to one vendor, your Purchasing Department will have fewer purchase orders and invoices. And a single point of contact for BP cuffs is created, providing you with a dedicated implementation team.

Products included in this catalog are not made with natural rubber latex.

## BP Cuff Adapters

### Connecting BP Cuff to DINACLICK and Neo-Snap Connectors



**Y, 1 Male Bayonet to DINACLICK™ Female Recept**  
 Compatibility: Philips/HP, Siemens, Datascope, Colin  
 Part # 2059291-001 (Qty: 10)



**Y, 1 Female Submin to DINACLICK Female Recept**  
 Compatibility: Protocol  
 Part # 2059291-002 (Qty: 10)



**Y, T Male Submin to DINACLICK Female Recept**  
 Compatibility: Nihon Kohden  
 Part # 2059291-003 (Qty: 10)



**Y, 1 Unfitted to DINACLICK Female Recept**  
 Compatibility: Unfitted  
 Part # 2059291-004 (Qty: 10)



**X, 2 Male Bayonet to DINACLICK Female Recept**  
 Compatibility: User Modified  
 Part # 2059292-001 (Qty: 10)



**X, 2 Female Submin to DINACLICK Female Recept**  
 Compatibility: GE Marquette  
 Part # 2059292-002 (Qty: 10)



**X, 2 Mated Submin to DINACLICK Female Recept**  
 Compatibility: GE Datex-Ohmeda  
 Part # 2059292-003 (Qty: 10)



**X, 2 Male Screw to DINACLICK Female Recept**  
 Compatibility: DINAMAP  
 Part # 2059292-004 (Qty: 10)



**Y, 1 Male Slip Luer to Neo-Snap\* Male Plug**  
 Compatibility: Philips/HP, Siemens, Datascope, Spacelabs  
 Part # 2059290-001 (Qty: 10)



**Male Micro Plug to 2 Neo-Snap Male Plug**  
 Compatibility: Spacelabs  
 Part # 2059290-002 (Qty: 10)

# BP Cuff Adapters

## Connecting BP Cuff to DINACLICK and Neo-Snap Connectors



**Male Micro Plug to 1 Neo-Snap Male Plug**

Compatibility: Spacelabs

Part # 2059297-001 (Qty: 10)



**Y, 1 Male Barb 1/8" to Neo-Snap Male Plug**

Compatibility: Philips, Welch Allyn

Part # 2059290-003 (Qty: 10)



**Y, 1 Male Barb 3/32" to Neo-Snap Male Plug**

Compatibility: Philips

Part # 2059290-004 (Qty: 10)

## Connecting BP Cuff to Luer



**Y, 1 Male Slip Luer to 2 Female Slip Luers-Neonatal**

Compatibility: Philips/HP, Siemens, Datascope, Spacelabs


Part # 330049 (Qty: 1)



**1 Male Micro Plug to 1 Female Locking Luer-Neonatal**

Compatibility: Spacelabs

Part # 330088 (Qty: 10)



**Y, 1 Male Micro Plug to 2 Female Locking Luers-Neonatal**

Compatibility: Spacelabs

Part # 330078 (Qty: 10)

## Connecting BP Cuff to 2 Female Screw Connectors



**Y, 1 Male Bayonet to 2 Female Screw Connectors**

Compatibility: Philips/HP, Siemens, Datascope, Colin

Part # 330046 (Qty: 1)



**Y, 1 Unfitted to 2 Female Screw Connectors**

Compatibility: Unfitted

Part # 330048 (Qty: 1)

## Connecting BP Cuff to 2 Unfitted



**Y, 1 Unfitted to 2 Unfitted**

Compatibility: Unfitted

Part # 2104081-001 (Qty: 1)

## Connecting BP Cuff to 2 Male Submins



**Y, 1 Unfitted to 2 Male Submins**

Compatibility: Unfitted

Part # 330084 (Qty: 10)



**Y, 1 Male Bayonet to 2 Male Submins**

Compatibility: Philips/HP, Siemens, Datascope, Spacelabs

Part # 330080 (Qty: 10)




**X, Mated Pair Submins to 2 Male Submins**

Compatibility: GE Datex-Ohmeda

Part # 330085 (Qty: 10)


# BP Cuff Connectors



**5/32" ID Barbed to  
1/8" ID Barbed, Plastic**  
Compatibility: Airhose Connector  
**Part # 300671 (Qty: 10)**



**5/32" ID Barbed  
to 5/32" ID Barbed, Plastic**  
Compatibility: Airhose Connector  
**Part # 300670 (Qty: 10)**



**Female Submin to  
5/32" ID Barbed, Plastic**  
Compatibility: GE, GE Marquette,  
GE Datex-Ohmeda  
**Part # 330092 (Qty: 10)**



**Male Screw to  
5/32" ID Barbed, Plastic**  
Compatibility: GE DINAMAP  
**Part # 300664 (Qty: 10)**



**5/32" ID Barbed to  
Female Bayonet, Metal**  
Compatibility: Philips/HP, Siemens,  
Datascope, Colin  
**Part # 330060 (Qty: 10)**



**Male Bayonet to  
3/16" ID Barbed, Metal**  
Compatibility: Philips/HP, Siemens,  
Datascope, Colin  
**Part # 330059 (Qty: 10)**



**2T DINACLICK Female Recept to  
1/8" ID tubing, dual lumen hoses**  
Compatibility: DINACLICK Male Plug  
**Part # 2059294-001 (Qty: 10)**



**1T DINACLICK Female Recept for  
single lumen hoses**  
Compatibility: DINACLICK Male Plug  
**Part # 2055630-001 1/8" tube  
(Qty: 10)**



**1T DINACLICK Female Recept for  
single lumen hoses**  
Compatibility: DINACLICK Male Plug  
**Part # 2055630-002 5/32" tube  
(Qty: 10)**



**5/32" ID Barbed to Female Slip Luer**  
Compatibility: Universal Neonatal  
**Part # 330068 (Qty: 10)**



**2T Neo-Snap Male Plug**  
Compatibility: Neo-Snap Female  
**Part # 2059293-001 (Qty: 10)**

## BP Cuff Connectors



**Male Submin to  
1/8" ID Barbed, Plastic**  
Compatibility: GE, GE Marquette,  
GE Datex-Ohmeda  
**Part # 330091 (Qty: 10)**



**Male Submin to  
5/32" ID Barbed, Plastic**  
Compatibility: GE, GE Marquette,  
GE Datex-Ohmeda  
**Part # 330090 (Qty: 10)**



**Mated Pair Submin to  
1/8" ID Barbed, Plastic**  
Compatibility: GE Datex-Ohmeda  
**Part # 330058 (Qty: 10)**



**Female Submin to  
1/8" ID Barbed, Plastic**  
Compatibility: GE,  
GE Marquette, Protocol  
**Part # 330064 (Qty: 10)**



**1/8" ID Barbed to  
Female Screw, Metal**  
Compatibility: GE DINAMAP  
**Part # 300619 (Qty: 2)**

## Wash Plugs and O-Rings

GE Healthcare provides wash plugs that block water from entering the BP cuff during cleaning.



**Washing Plugs for  
Submin Connector Cuffs**  
Compatibility: GE,  
GE Marquette, Protocol  
**Part # 330072 (Qty: 50)**



**Washing Plugs for  
Screw Connector Cuffs**  
Compatibility: GE DINAMAP  
**Part # 300877 (Qty: 50)**



**O-Ring Replacement Set for  
Air Hose Instrument Connector**  
Compatibility: (330073) and Air Hoses  
(107365, 107368, 107366)  
**Part # 330077 (Qty: 20)**



**Washing Plugs for Neo-Snap  
Female Recept and DINAClick  
Male Plug**  
Compatibility: GE  
**Part # 2059295-001 (Qty: 50)**

## For More Information

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at [cuffs.gehealthcare.com](http://cuffs.gehealthcare.com)

## Imagination at work

© 2018 General Electric Company - All rights reserved.

General Electric reserves the right to make changes in specifications and features, or discontinue the product or service described at any time, without notice or obligation. This does not constitute a representation or warranty or documentation regarding the product or service featured. Illustrations are provided for informational purposes, and your configuration may differ.

GE and GE Monogram, CRITIKON, DINACLICK are trademarks of General Electric Company.

[www.gehealthcare.com](http://www.gehealthcare.com)

Colin is a registered trademark of Colin Medical Instruments.

Datascope is a registered trademark of Datascope Corp.

HP is a registered trademark of Agilent Technologies.

Nihon Kohden is a registered trademark of Nihon Kohden Corporation.

Philips is a registered trademark of Royal Philips Electronics.

Protocol is a registered trademark of Welch Allyn, Inc.

Siemens is a registered trademark of Siemens Medical Solutions.

Spacelabs is a registered trademark of Spacelabs Medical Systems.

GE Healthcare, a division of General Electric Company.

GE Medical Systems *Information Technologies*, Inc., doing business as GE Healthcare.

