



GE Healthcare

Clinical Accessories that support **Infection Prevention & Control**

Infection Prevention and Control Policy is to ensure that the Hospital's provide an environment and system of care which minimises the risk of infection to patients, staff and visitors.

GE Accessories offer a reliable disposable solution that delivers value.

- Single patient use accessories **reduce the risk of cross-contamination** and help **lower the risk of costly equipment associated infections**
- No sterilization or assembly required, which could lead to **time and cost savings**
- **Help reduce the risk of nosocomial infections**

INTENSIVE CARE

Single-Patient-Use Valve



Flow Sensor Only
1505-3231-001



Housing only
(includes water trap)
1505-8568-001



Whole assembly
(Housing and Flow sensor)
1505-3848-001



Inspiratory Safety Guard
2083208-001 (10 pcs/Pkg)

NIBP cuffs



Yellow Critikon™ Classic-Cuf*



Soft-Cuf*



Radial-Cuf*

Disposable SpO₂ Sensor



TruSignal AllFit Sensor
TS-AF-10 (10 pcs/Pkg)
TS-AF-25 (25 pcs/Pkg)



TruSignal Adult/Pediatric Sensor
TS-AP-10 (10 pcs/Pkg)
TS-AP-25 (25 pcs/Pkg)



D-fend+
10/case
881319-HEL¹



D-fend Pro+
10/case
M1200227¹



Gas Sampling line
2m - 10/case
2097307-001

Airway Gas Sampling Set, Disposable Adult/Pediatric, TPE, 2m/10ft
20 pcs/Pkg
8090121521

D-lite+ sensor
50/case
896952

PERIOPERATIVE CARE

Gas Sampling line & Flow sensor



Gas Sampling line
2m - 10/case
2097307-001



D-lite+ sensor
50/case
896952



D-fend
10/case
876446-HEL²



D-fend Pro
10/case
M1182629²



Adult HEPA pleated paper hydrophobic bacterial viral filter and Luer port
2106570-008

NIBP cuffs

| SOFT-CUF | | | GE Healthcare | | | | Non-GE |
|-------------|------------|--------------------|------------------------|----------------------------|-------------------------|-------------------------------|--------------------------|
| Size | Color | Limb Circumference | 2-Tube Screw Connector | 2-Tube DINACLICK Connector | 2-Tube Submin Connector | 2-Tube Mated Submin Connector | 1-Tube Bayonet Connector |
| Small Adult | Light Blue | 17-25 cm | 2502 ³ | SFT-A1-2A ³ | 2403 ³ | 2452 | 2362 ³ |
| Adult | Navy | 23-33 cm | 2503 ³ | SFT-A2-2A ³ | 2404 ³ | 2453 | 2363 ³ |
| Adult Long | Navy | 23-33 cm | 2604 ³ | SFT-A2-2A-L ³ | 2116 ³ | 2454 | 2140 ³ |
| Large Adult | Rose | 31-40 cm | 2504 ³ | SFT-A3-2A ³ | 2405 ³ | 2455 | 2364 ³ |

| RADIAL-CUF | | GE Healthcare | | | |
|------------|-----------------------|------------------------|----------------------------|-------------------------|-------------------------------|
| Color | Forearm Circumference | 2-Tube Screw Connector | 2-Tube DINACLICK Connector | 2-Tube Submin Connector | 2-Tube Mated Submin Connector |
| Purple | 26-36 cm ⁴ | SFT-F1-2E | SFT-F1-2A | SFT-F1-2C | SFT-F1-2D |

Note: GE Critikon RADIAL-CUF BP cuffs come 20 in a box.

| YELLOW CLASSIC-CUF | | | GE Healthcare | | |
|--------------------|-------------|--------------------|------------------------|----------------------------|-------------------------|
| Size | Color | Limb Circumference | 2-Tube Screw Connector | 2-Tube DINACLICK Connector | 2-Tube Submin Connector |
| Small Adult | Yellow/Blue | 17-25 cm | 2607 | CLA-A1-2A-X | 2672 |
| Adult | Yellow/Blue | 23-33 cm | 2602 | CLA-A2-2A-X | 2673 |
| Large Adult | Yellow/Blue | 31-40 cm | 2642 | CLA-A3-2A-X | 2674 |

GE masks and breathing circuits are also available at your service, please contact your local sales representatives for further information.

gehealthcare.com

© 2020 General Electric Company – All rights reserved.

GE Healthcare reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Healthcare representative for the most current information.

GE, GE Monogram are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company.

¹ Recommended to change water trap for each new patient, whenever the host device prompts you to do so, or every 24 hours.

² Change whenever the host device prompts you to do so. The maximum lifetime of a water trap is two months.

³ Has a respective part number with an "-INT" suffix, denoting that the Instructions for Use (IFU) is available in multiple languages.

⁴ Forearm measurement